



Test and track hydration anytime, anywhere

Hydration Testing System



	Item Code: HTS-01 Standard HTS <i>Recommended for 1-15 individuals.</i>	Item Code: HTS-01E Travel HTS <i>Recommended for 1-15 individuals.</i>	Item Code: HTS-02 Gold Package <i>Recommended for 30-60 individuals.</i>	Item Code: HTS-03 Platinum Package <i>Recommended for 60-100 individuals.</i>
MX3 Lab Pro	1	1	3	5
Hydration Test Strips	120 (3 boxes of 40)	120	360	600
Additional Test Strips	-	-	600	1000
Dispenser	1	1	3	5
USBC Charging Cable	1	1	3	5
Travel Case	-	1	3	5
Wrist Straps	-	-	3	5
MX3 Mobile App	Included	Included	Included	Included
MX3 Web Portal	Included	Included	Included	Included
Price (USD)	A\$ 1360 + GST	A\$ 1360 + GST	A\$ 4080 + GST	A\$ 6800 + GST
Bundle Savings	-	-	A\$ 416 + GST	A\$ 694 + GST

Hydration Test Strips

Accessories



Item Code: SH-01
Single (40 strips)
A\$ 32 + GST
(A\$ 0.80/test)

Item Code: SH-02
5-Pack (200 strips)
A\$ 150 + GST
(A\$0.75/test)

Item Code: SH-03
50-Pack (2000 strips)
A\$ 1360 + GST
(A\$0.68/test)

Test Strips Dispenser
A\$ 36 + GST

Travel Case
A\$36 + GST



USB-C Charging Cable
AUD 18 + GST



Wrist Strap
AUD 6.80 + GST

*All prices excludes GST and are subject to change.

Each Hydration Testing System Kit Includes

MX3 Lab Pro



This lightweight, handheld device can go anywhere. Simply insert a test strip and take a measurement for immediate hydration status. Data automatically pushes to the App and into the cloud for Dashboard access, if desired.

Hydration Test Strips



These disposable test strips need just a micro-liter of saliva from the tip of the tongue. Results are delivered in seconds. Test strips can be refilled monthly or as needed. (qty 120)

Dispenser



Lightweight and portable, the dispenser stores strips and provides a hygienic, hands-free way to load the LAB Pro.

MX3 App



The App is available for a free download on Apple and Android devices. (Included)

MX3 Online Dashboard



Access the Dashboard via the online portal to review user profiles and all measurement data. The Dashboard provides monitoring and tracking over time, with the analytics you need to help personalize hydration strategies and identify organization-wide trends. (Included)